

Motion of a Thrown Ball

We can measure the distance of a ball thrown vertically, analyze its speed, average acceleration, potential energy, and kinetic energy.

Fundamental Concept

When a ball is thrown straight up, its speed decreases until it reaches its peak height. At this point, the kinetic energy becomes zero while the potential energy increases to its maximum. After reaching the highest point, the ball falls back down under the influence of gravity, increasing in speed.



| Height | Speed | Potential energy | Kinetic energy | Acceleration |
|----------|----------|------------------|----------------|--------------|
| Increase | Decrease | Increase | Decrease | g |







| Height | Speed | Potential energy | Kinetic energy | Acceleration |
|----------|----------|------------------|----------------|--------------|
| Decrease | Increase | Decrease | Increase | g |


Experiment

Materials Needed

Interface, Science# program, Smart device, Motion sensor, Basketball

Interface Setup

1.  Run the Science# program.
2. Connect the motion sensor to the interface.
3. Press  to set up the experimental environment as shown below or press   for automatic setup.



Experiment Setting

Data collection method

☒ Auto collection

☐ Manual collection

☐ data collect as absolute value

Chart type

☒ Line chart

☐ Bar chart

☐ X-Y chart

Data on the X-axis :


Data collecting interval


Experiment by time


Data count: 3000

☐ Display the current time on the x-axis

Data Collection

1. Place the motion sensor on the ground facing the ceiling.
2. Stand holding the basketball above the sensor.
3. Press  to start data collection.
4. Throw the basketball straight up. Ensure the ball stays within the range of 0.15–2m from the motion sensor.



5. When the motion of the basketball stops, press  to end data collection.

Data Analysis

Recording Data

1. Draw a time-distance graph of the ball thrown vertically and explain the distance changes according to the ball's movement direction.

| Motion of the Ball | Going Up | Coming Down |
|--------------------|----------|-------------|
| Distance Change | | |

Data Application

1. Use [Analyze] - [Differentiate] to draw a time-speed graph of the ball thrown vertically and explain the speed changes according to the ball's movement direction.

| Motion of the Ball | Going Up | Coming Down |
|--------------------|----------|-------------|
| Speed Change | | |

2. Analyze the time-speed graph of the ball thrown vertically using [Analyze] - [Linear $f(x)=Ax+B$] to draw the graph and find the average acceleration (A).

| Average Acceleration of the Vertically Thrown Ball (m/s ²) |
|--|
| |

3. Analyze the time-distance graph of the ball thrown vertically using [Analyze] - [Conservation of Mechanical Energy] - [Potential Energy $E_p=mgh$] to draw the potential energy graph, and use [Analyze] - [Conservation of Mechanical Energy] - [Kinetic Energy $E_k=(1/2)mv^2$] to draw the kinetic energy graph.

| Motion of the Ball | Going Up | Coming Down |
|-------------------------|----------|-------------|
| Potential Energy Change | | |
| Kinetic Energy Change | | |

